



## Thermo Scientific™ TruNarc™ Delta Handheld Narcotics Analyzer

Rapid, contactless, and reliable identification of narcotics for presumptive testing in an easy-to-use field device.

The Thermo Scientific™ TruNarc™ Delta Handheld Narcotics Analyzer expands on the legacy of the Thermo Scientific™ TruNarc™ Handheld Narcotics Analyzer, which has been trusted for over a decade by law enforcement and government agencies worldwide.

- **Comprehensive Testing:** Tests for over 1,200 substances, including narcotics, precursors, and cutting agents.
- **Continuously Updated Library:** Regular updates to include emerging narcotics threats.
- **Contactless Operation:** Requires no direct contact with most substances, enhancing safety.
- **Performance Enhancements for Faster Results:** Provides clear, real-time results for presumptive evidence.
- **Advanced Type H2 Kit:** Capable of testing fluorescent substances and low-concentration mixtures.
- **Automated Record Keeping:** Generates records with complete scan results, including time-and-date stamps, to expedite prosecution.

### Specifications

Weight	1.27 lbs (0.575 kg)
Size	6.5 in x 4.1 in x 2 in (16.5 cm x 10.4 cm x 5.1 cm)
Library	Controlled substances, cutting agents, and precursors
Data export formats	CSV, SPC, PDF, SCZ
Battery	Rechargeable internal 3.6V battery (10hrs battery life), user replaceable. DC wall adapter, 5V@3A output USB-C
Operating temperature	-10°C to +50°C (14°F to 122°F)
Language configurations	English
Computer administration	TruNarc Admin software connected via USB-C
Reachback support	Automated spectral analysis with additional feedback from staff scientists as needed

■ To learn more about the TruNarc or schedule a demo, please visit [safewareinc.com](https://safewareinc.com)



800.331.6707 | [info@safewareinc.com](mailto:info@safewareinc.com)

thermo scientific